NuBODY

NÜFACE®

USER MANUAL

ENSURE YOUR BEST RESULTS!

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Zumbido energizante (no se muestra)

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INTENDED USE

The NuBODY Skin Toning Device is intended for body skin stimulation and is indicated for over-the-counter cosmetic use.

IMPORTANT INFORMATION

Use NuBODY only as described within this User Manual. The risks and dangers of using NuBODY in any way that does not follow this User Manual might be unknown.

Your NuFACE NuBODY Device has a unique serial number located on the back of the device. Record your serial number and register your device online at www.MyNuFACE.com. A registration link can be found at the bottom of the webpage.

CONTRAINDICATIONS FOR USE

- Do not use this device if you have a history of epilepsy or seizures.
- Do not use if you have had recent surgery or other surgical procedure – consult your healthcare provider.
- Do not use this device if you have a heart pacemaker, heart rhythm stimulator placed in your body (implanted defibrillator), or other electronic device put inside your body. Such use could cause electric shock, burns, electrical interference, or death.
- Do not use if you have any active cancer or over any suspicious or cancerous lesion – consult your healthcare provider.
- Do not use this device if you are pregnant because the safety of electrical stimulation during pregnancy has not been established.
- Do not use on infants and children under the age of 18, and keep device out of the reach of children.
- This device is not intended for use by persons (including children) with reduced physical, sensory, or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the device.
- Consult your healthcare provider if any of the contraindications above apply to you.

⚠ WARNINGS

- If you are under the care of a physician, consult your physician before using this device.
- Do not operate NuBODY device in the presence of electronic monitoring equipment (e.g., heart monitors, heart tracing) or in close proximity (e.g., 1 m) to shortwave, microwave or electromagnetic therapy equipment as it may cause the medical equipment to malfunction. Such equipment may produce instability in NuBODY device output and cause burns at the stimulation site.
- Stimulation should not be applied over or near cancerous wounds or open wounds such as rashes, red, swollen, infected, or inflamed areas or skin eruptions, e.g., phlebitis, thrombophlebitis, broken capillaries, varicose veins, etc.
- Stimulation should not be applied over your neck as this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Apply stimulation only to normal, intact, clean, healthy skin.

Consult your healthcare provider:

- if you experience an adverse reaction from using the device and discontinue use.
- before using the device if you have painful areas.

Do not apply stimulation / use the device:

- across your chest (thorax area) may increase the risk of cardiac fibrillation or other rhythmic disturbances, which could be lethal.
- on the head, face (including directly on the eyes and covering the mouth), neck (especially the carotid sinus), groin, upper back, chest area or crossing over the heart.
- over areas of skin that lack normal sensation.
- while sleeping, bathing or showering.
- when there is a tendency to bleed following acute injury or breaks in bones.

⚠ WARNINGS (CONTINUED)

- on the head, face (including directly on the eyes and covering the mouth), neck (especially the carotid sinus), groin, upper back, chest area or crossing over the heart.
- while driving, operating machinery, or during any activity in which uncontrollable muscle twitching may put the user or the public at unnecessary risk of injury.
- if it is not functioning normally, or not operating properly due to damage caused by dropping or accidental submersion in water.
- if you have a body-worn medical device such as an insulin pump.

Discontinue use:

- If the device begins to overheat, turn the device OFF immediately.
- Keep all electrical appliances (including the NuBODY device) away from water (including baths, showers, and sinks).
- Never drop or insert any object into any openings on the device.
- If you suspect the device to be malfunctioning, or not operating properly please contact NuFACE.
- Use of accessories, transducers and cables other than those specified, except for transducers and cables sold by the manufacturer of the NuBODY device as replacement parts for internal components, may result in increased emissions or decreased immunity of the system.
- Keep this device out of the reach of children and pets.
- Use of controls or adjustments, or the performance of procedures other than those specified herein, may result in injury.
- Do not disassemble the device as this may cause damage, malfunction, electrical shock, fire, or injury. There are no user-serviceable parts inside. **Note: Disassembly voids the warranty.**

PRECAUTIONS

- The NuBODY device is designed for cosmetic use only and for individuals in good health. If this is not the case for you, please consult your doctor before use.
- The long-term effects of transcutaneous electrical stimulators are unknown.
 - If you have suspected or diagnosed heart disease, you should follow precautions recommended by your physician.
 - Use caution if stimulation is applied over the menstruating uterus.
 - You may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (e.g., the NuFACE Microcurrent Activator). Irritation can usually be reduced by using an alternate conductive medium. If you suspect that you have had a negative reaction from using the electrical stimulation device, consult your doctor immediately.
 - Skin irritation and burns beneath the electrodes have been reported with the use of other powered muscle stimulators. If you experience skin irritation or burning, discontinue use immediately.
 - If using NuBODY causes persistent pain or discomfort, discontinue use immediately.
 - A slight tingling sensation may occur while using NuBODY. If you feel the current, apply more NuFACE Microcurrent Activator. Lowering the intensity may reduce or eliminate this sensation.
 - The lithium-ion battery of NuBODY is not removable by the user. If the device is to be disposed of, do not put it into a fire. Dispose of the device per local regulations.
 - Inspect the device and do not use the device if you see any signs of internal battery leaks.
 - Do not store NuBODY in the sun or on a hot surface. High temperatures may be damaging.
 - The device should be operated, transported and stored in the manner specified in this manual.
 - Failure to use and maintain the device in accordance with the instructions outlined in this manual will void the warranty.
 - Only use the external power supply provided with the NuBODY with the following specifications: **5VAC 1A**. The power adapter is marked NuFACE. Always position the power adapter so it is accessible in case it is needed to be disconnected from the mains.

NuBODY is intended for cosmetic use on the body as described. Any consequence resulting from misuse, application to other parts of the body or face, connection to improper voltage sources, dirty conductive solution or spheres or any other improper application is neither the responsibility of NuFACE nor its affiliates.

The NuFACE NuBODY has been designed to provide freedom from unacceptable risk when used under normal conditions and reasonably foreseeable hazards. Even in the event of complete loss of treatment function, NuBODY does not pose a risk to the User.

BEFORE USING YOUR NEW DEVICE

- 1) Remove your NuBODY and power adapter from the box.
- 2) Place your NuBODY on a counter or table top. Keep away from any source of water.
- 3) Connect the power adapter to the NuBODY. Plug adapter into the wall outlet to charge your device. **Note: NuBODY cannot be used while plugged into the outlet and will not turn ON while charging.**
- 4) **Important: Charge your NuBODY for approximately 16 hours before the first use.** The blue LED lights will blink, showing that your NuBODY is charging. All 3 blue LED lights will stay lit when NuBODY is fully charged.
- 5) **Important: Test your skin to make sure you are not sensitive to the NuFACE Microcurrent Activator.** To do so, apply a dime-sized amount of the NuFACE Microcurrent Activator on your forearm. Leave for 15–30 minutes before washing it off. If an allergic reaction (rash, redness, etc.) develops, immediately contact NuFACE for recommendations on alternate NuFACE Microcurrent Activators. You should not initiate treatment with the NuBODY device until your reaction to the NuFACE Microcurrent Activator has been resolved.
- 6) When NuBODY is ON it will emit an energizing hum. The default intensity is set at the minimum level.
- 7) Always ensure that all four [4] of the spheres remain in direct contact with the treatment area. Always glide upward with medium pressure without pulling or dragging the skin.
- 8) NuBODY emits a tonal beep every five (5) seconds to indicate when to move to the next position if desired. This beep will change to two (2) tonal beeps when the battery needs recharging. **Note: NuBODY will NOT contract the muscles. The current is soft, soothing, and mild when used correctly.**

DEVICE FEEDBACK — INDICATES

Ascending Beeps	Powered ON
Descending Beeps	Powered OFF
Tonal Beep	Intensity Adjusted
Dual-Tone Beep	End of Glide Interval

Two series of tonal beeps when device powers OFF indicates the battery requires charging

NUFACE MICROCURRENT ACTIVATOR

The essential first step to your microcurrent treatment to ensure optimal results. **Important: Always apply Microcurrent Activator when using the NuBODY Device.** Refer to packaging for use by date.

TREATMENT AREAS

The NuBODY Device is designed with four treatment spheres designed to stimulate the larger, denser areas of the body. Follow instructions below for best results and refer to your Quick Start for more detailed images.

TREATMENT TECHNIQUES

GLIDES:

To glide, place NuBODY at the starting point of a treatment area. When it beeps, slowly glide the device in an upward motion using medium pressure to the endpoint until the beep sounds (5 seconds).

HOLDS (Recommended for stubborn areas):

To hold, place NuBODY on the targeted area and hold in place for at least 2-3 beeps (10-15 seconds). Do not glide. You may treat multiple stubborn areas within the treatment area for up to 5 minutes.

HOW TO USE YOUR NUBODY DEVICE

Get ready to Smooth and Tone.

The instructions below are designed to provide the most effective results when followed correctly. **Note: For best results, you may choose to focus on a specific area of concern within the treatment area for the full 5 minute treatment cycle and then repeat over other areas as necessary.**

- 1) Cleanse the selected treatment area with an oil-free cleanser. **Note: Any cream, lotion, self tanner, and/or oil-based products on the skin can block the microcurrent.** **Note: Shave any excessive hair in the treatment area as it can create a tingling sensation.**
- 2) **Important: Apply a mask-like layer of NuFACE Microcurrent Activator in sections within the selected treatment area.** **Note: Always apply a NuFACE Microcurrent Activator when using your NuBODY Device. Reapply if skin becomes dry during treatment to avoid any uncomfortable sensations.**
- 3) Unplug the power adapter from your NuBODY device.
- 4) Turn ON NuBODY by pressing the ON/OFF button. Ascending tonal beeps will sound. Select your preferred intensity level.
- 5) Adjust the intensity level up or down by pressing the power button until you reach your comfort level. As you move between the three (3) intensity levels, NuBODY will emit a tonal beep and the number of blue LED lights will change.

Intensity Levels:

- Minimum Intensity: First Light Illuminated
- Medium Intensity: First 8 Second Lights Illuminated
- Maximum Intensity: All Three Lights Illuminated

- 6) Begin treatment by gliding the device in an upward motion as shown in the illustration below.
- 7) For stubborn areas, such as dimpled or uneven skin, hold the device on the target area for 2-3 beeps to provide a more targeted treatment.
- 8) Treat for 5 minutes until the device automatically powers OFF. You may continue to treat other areas for an additional 5 minutes per area.
- 9) After treatment is complete, massage any excess Microcurrent Activator into skin for additional skincare benefits.
- 10) Turn OFF the NuBODY Device after treatment. A descending audible beep is emitted notifying you that your device is turning OFF. See the "CLEANING" section for proper cleaning instructions. **Note: The NuBODY Device will automatically shut off after 5 minutes of use.**

Important: Download the NuFACE APP for additional how-to instructions, and customized treatment options. You can download it from the App store or the Google Play Store. To pair your device, in the App, click on 'Add Device' and follow the on-screen instructions to connect. The NuFACE app is available in many countries around the world and supports Chinese, Spanish, German, French, and Korean languages.

TREATMENT FREQUENCY

- Use 5 minutes a day per treatment area.
- For the first 60 days, use 5 times per week.
- After 60 days, use 2-3 times per week for maintenance.
- Allow 24 hours between treatments.

RESULTS

Results may vary based on skin type and skin condition. Best results are achieved with regular use. **This is not a weight loss product.**

CARING FOR YOUR DEVICE

Cleaning: Slightly dampen a soft, lint-free cloth with water or rubbing alcohol to clean your device after every treatment. **Do not immerse device in water at any time.**

Charging: Before first use, charge your NuBODY for 16 hours. For ideal battery life, the device should be continuously charged between treatments.

Battery Care: It is advised to store NuBODY connected to the plugged-in power adapter after each use. **Note: The rechargeable battery of NuBODY is not removable by the user.** Use caution prior to using the device if it has been stored or has not been used for extensive periods. Inspect the device and do not use the device if any signs of internal battery leaks are evident. **Important: If the device is to be disposed of, do not put into a fire. Dispose of the device and its batteries per local regulations.**

Troubleshooting: If your NuBODY is not turning ON, make sure the battery is fully charged. If not, charge your device for a full 16 hours.

If your device is still not working, if you need assistance, if you need to report unexpected events, or if you have any questions, please call NuFACE Client Services at 888-360-9171

WARRANTY

NuFACE® warrants that all NuFACE® NuBODY devices are free from manufacture defects in materials and workmanship for a period of one (1) year from the date of your original purchase. Our obligation under this Warranty is limited to replacement of the product which will not extend the period of your Warranty. This Warranty does not cover damage covered by abuse, improper care, application of harmful chemicals, product that may have been altered in any way, loss of the device and product purchased from an unauthorized seller. NO OTHER WARRANTY EXPRESSED OR IMPLIED, WHETHER OF FITNESS FOR A PARTICULAR PURPOSE OR OF MERCHANTABILITY OR OF ANY OTHER KIND, SHALL EXIST IN RESPECT TO SUCH PRODUCT.

For Devices Purchased in the US: For more information about your warranty or to initiate a warranty claim and receive a return authorization number, contact our Client Services Team via our website at www.MyNuFACE.com/warranty or by calling 888-360-9171. Must show receipt from an authorized seller as proof of purchase.

For Devices Purchased Outside of the US, International Returns and International Warranties: For all product purchased outside of the United States, warranty and returns are per your place of purchase—spa, retailer or distributor. Please consult your receipt or place of purchase for their return and warranty policy.

Conforms to IEC 60529, 60601-1, 60601-1-2, 60601-1-6, 60601-2-10, 62366

SYMBOLS

One or more of the following markings may appear on the NuBODY Device, User Manual or packaging.

	European Conformity		Comply with European Restriction of Hazardous Substances Directive 2002/95/EC
	European Representative		Product must be disposed of in accordance with the WEEE directive
	Part Number		FCC-approved equipment authorization mark
	Operator's Manual; Operating Instructions		Class II double insulated electrical device
	Refer to Instruction Manual/Booklet		

SYMBOLS (CONTINUED)

	Manufacturer
	Power On/Off
	Detachable power supply port
	Type BF applied part, which is equipment that provides a particular degree of protection against electric shock. Particularly regarding allowable leakage current and of the protective earth connection.

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All Rights Reserved. / Tous droits réservés. / Todos los derechos reservados. / Alle Rechte vorbehalten. / 모든 권리 보유.
NuFACE® and Trinity® are registered trademarks of the Carol Cole Company Dba NuFACE.
Medical Device Safety Service GmbH